

NA

InVivoMAb Anti-Influenza A virus Neuraminidase/NA (1G01)

Catalog No.	VVV24001A VVV24001B	Quantity:	100 µg 1.0 mg
Alternate Names:	Neuraminidase, NA		
Description:	Neutralizing human monoclonal antibody (clone 1G01) which targets Neuraminidase from Influenza A virus (A/Minnesota/11/2010(H3N2))		
UniProt ID (target):	C3W5S3		
Source:	Mammalian cells		
Specificity:	Influenza A virus (A/Minnesota/11/2010(H3N2)) Neuraminidase		
Isotype:	Human IgG1		
Clone:	1G01		
Purification:	Protein A or G affinity chromatography from cell culture supernatant		
Purity:	> 95% by SDS-PAGE and SEC-HPLC		
Formulation:	Sterile-filtered 0.01M PBS, pH 7.4		
Endotoxin:	< 0.01 EU/µg		
Concentration:	1.0 mg/ml, lot specific		
Applications:	Functional studies, Neutralization, ELISA		
Storage & Stability:	Store at 2-8°C for immediate use within 1-2 weeks or aliquot the stock solution and store at -20°C to -80°C up to 1 year. Avoid repeated freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences®
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com