

N

Recombinant HMPV N/Nucleoprotein, N-His Tag

| | | | |
|---------------------------------|---|------------------|------------------|
| Catalog No. | YVV30001A YVV30001B | Quantity: | 100 µg 1.0 mg |
| Alternate Names: | Nucleoprotein, Protein N, Nucleocapsid protein | | |
| Description: | Recombinant Human metapneumovirus (strain CAN97-83) (HMPV) nucleoprotein Met1-Glu394 with a N-terminal His tag expressed in E. coli. | | |
| UniProt ID: | Q6WBA1 | | |
| Protein Construction: | Met1-Glu394 for HMPV nucleoprotein with a N-terminal His tag | | |
| Source: | E. coli | | |
| Formulation: | Lyophilized from sterile-filtered PBS, pH 7.4, 0.02% NLS, 1 mM EDTA, 4% trehalose, 1% mannitol | | |
| Molecular Weight: | 45.70 kDa, predicted | | |
| Purity: | > 90 % as determined by SDS-PAGE. | | |
| Reconstitution: | Centrifuge vial prior to opening. Add sterile distilled water to a concentration of 1.0 mg/mL and gently pipette the solution up and down the sides of the vial. DO NOT VORTEX. Allow several minutes for complete reconstitution. | | |
| Storage & Stability: | Stable for up to 1 year from date of receipt at -20°C to -80°C. After reconstitution, stable for 1 week at 2-8°C, or store working aliquots at -20°C to -80°C. Avoid repeated freeze-thaw cycles. | | |

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences[®]
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com