

Recombinant WEEV Spike glycoprotein E1, N-His Tag

Catalog No.	YVV34402A YVV34402B	Quantity:	100 µg 1.0 mg
Alternate Names:	Spike glycoprotein E1, E1 envelope glycoprotein		
Description:	Recombinant Spike glycoprotein E1 (Ile810-Asn1205) from Western equine encephalitis virus (WEEV) with a N-terminal His tag, was expressed in E. coli.		
UniProt ID:	P13897		
Protein Construction:	Spike glycoprotein E1 (Ile810-Asn1205) from Western equine encephalitis virus (WEEV) fused with a N-terminal His tag		
Source:	E. coli		
Formulation:	Lyophilized from sterile-filtered PBS, pH 7.4, 0.02% NLS, 1 mM EDTA, 4% trehalose, 1% mannitol		
Molecular Weight:	45.08 kDa, predicted		
Purity:	> 90 % as determined by SDS-PAGE.		
Applications:	SDS-PAGE, Western Blot, ELISA, Immunogen		
Reconstitution:	Centrifuge vial prior to opening. Add sterile distilled water to a concentration of 1.0 mg/mL and gently pipette the solution up and down the sides of the vial. DO NOT VORTEX. Allow several minutes for complete reconstitution.		
Storage & Stability:	Stable for up to 1 year from date of receipt at -20°C to -80°C. After reconstitution, stable for 1 week at 2-8°C, or store working aliquots at -20°C to -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

