

## RHDV

### Recombinant RHDV VP60/Capsid protein VP60, C-His Tag

<b>Catalog No.</b>	YVV35601A YVV35601B	<b>Quantity:</b>	100 µg 1.0 mg
<b>Alternate Names:</b>	Genome polyprotein, Capsid protein VP60		
<b>Description:</b>	Recombinant Capsid protein VP60 (Thr2009-Val2344) from Rabbit hemorrhagic disease virus (RHDV) was fused with a C-terminal His tag and expressed in E. coli.		
<b>UniProt ID:</b>	Q3HNQ2		
<b>Protein Construction:</b>	Capsid protein VP60 (Thr2009-Val2344) from Rabbit hemorrhagic disease virus (RHDV) was fused with a C-terminal His tag		
<b>Source:</b>	E. coli		
<b>Formulation:</b>	Lyophilized from sterile-filtered PBS, pH 7.4, 0.02% SKL, 1 mM EDTA, 5% trehalose, 2% mannitol		
<b>Molecular Weight:</b>	35.91 kDa, predicted		
<b>Purity:</b>	> 90 % as determined by SDS-PAGE.		
<b>Applications:</b>	SDS-PAGE, Western Blot, ELISA, Immunogen		
<b>Reconstitution:</b>	<b>Centrifuge vial prior to opening.</b> Add sterile distilled water to a concentration of 1.0 mg/mL and gently pipette the solution up and down the sides of the vial. <b>DO NOT VORTEX.</b> Allow several minutes for complete reconstitution.		
<b>Storage &amp; Stability:</b>	Stable for up to 1 year from date of receipt at -20°C to -80°C. After reconstitution, stable for 1 week at 2-8°C, or store working aliquots at -20°C to -80°C. <b>Avoid repeated freeze-thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



**Cell Sciences**<sup>®</sup>  
65 Parker Street  
Unit 11  
Newburyport, MA 01950

Toll Free: 888-769-1246  
Phone: 978-572-1070  
Fax: 978-992-0298

E-mail: [info@cellsciences.com](mailto:info@cellsciences.com)  
Website: [www.cellsciences.com](http://www.cellsciences.com)