

emm3

## Recombinant *Streptococcus pyogenes* M protein, N-His Tag

<b>Catalog No.</b>	YXX21704A YXX21704B	<b>Quantity:</b>	100 µg 1.0 mg
<b>Alternate Names:</b>	M protein		
<b>Description:</b>	Recombinant <i>Streptococcus pyogenes</i> M protein, Leu63-Asn112 expressed in <i>E. coli</i>		
<b>UniProt ID:</b>	W0T370		
<b>Source:</b>	Expressed in <i>E. coli</i> , cloned from <i>Streptococcus pyogenes</i> protein M		
<b>Molecular Weight:</b>	14.35 kDa, Leu63-Asn112		
<b>Formulation:</b>	Lyophilized from PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% trehalose, 1% mannitol		
<b>Purity:</b>	> 90 % as determined by SDS-PAGE.		
<b>Tag:</b>	N-terminal His		
<b>Reconstitution:</b>	Reconstitute with sterile deionized water to prepare a stock solution of 0.1-1.0 mg/ml.		
<b>Storage &amp; Stability:</b>	Upon receipt, store lyophilized protein at -20°C to -80°C for up to 1 year. <b>Avoid repeated freeze-thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



**Cell Sciences**<sup>®</sup>  
65 Parker Street  
Unit 11  
Newburyport, MA 01950

Toll Free: 888-769-1246  
Phone: 978-572-1070  
Fax: 978-992-0298

E-mail: [info@cellsciences.com](mailto:info@cellsciences.com)  
Website: [www.cellsciences.com](http://www.cellsciences.com)