

emm

Recombinant European house dust mite Der p 10/ Tropomyosin Protein, N-His Tag

Catalog No.	YXX53301A YXX53301B	Quantity:	100 µg 1.0 mg
Alternate Names:	Tropomyosin, Der p 10		
Description:	Recombinant <i>Dermatophagoides pteronyssinus</i> (European house dust mite) Tropomyosin protein, [Met1-Tyr284] expressed in <i>E. coli</i>		
UniProt ID:	O18416		
Source:	<i>E. coli</i>		
Molecular Weight:	35.06 kDa		
Formulation:	Lyophilized from PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% trehalose, 1% mannitol		
Purity:	> 90 % as determined by SDS-PAGE.		
Tag:	N-terminal His		
Reconstitution:	Reconstitute with sterile deionized water to prepare a stock solution of 0.1-1.0 mg/ml.		
Storage & Stability:	Upon receipt, store lyophilized protein at -20°C to -80°C for up to 1 year. Avoid repeated freeze-thaw cycles.		



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.